



# Differential Electrochemical Mass Spectrometry (DEMS)


## SOLUTIONS FOR DISSOLVED GAS ANALYSIS AND OFF-GAS ANALYSIS IN ELECTROCHEMISTRY

### Key Features

- ▶ Differential Electrochemical Mass Spectrometry (DEMS) is an analytical technique that combines electrochemical half-cell experimentation with mass spectrometry
- ▶ It allows in situ mass resolved determination of gaseous or volatile electrochemical reactants, reaction intermediates and products in real time
- ▶ Hiden Analytical offer a range of DEMS cells with electrolyte/nanoporous sampling interface to the Hiden HPR-40 DSA Mass Spectrometer



Hiden Analytical DEMS Cell



▶ Vitreous Carbon – User coatable working electrode


▶ 4 Ports for additional electrodes

▶ Replaceable nanoporous membrane


▶ Interface to the Hiden HPR-40 DSA

## Novel DEMS Cell

Hidden Analytical have a new license agreement with The University of California, Berkeley Lab. Under the terms of the agreement Hidden Analytical will commercialize a novel differential electrochemical mass spectrometry (DEMS) cell for integration with the Hidden Analytical range of membrane inlet mass spectrometer systems. The new cell, developed by Ezra L. Clark and Prof. Alexis T. Bell, coupled with Hidden's differentially pumped mass spectrometer system provides for in-situ mass resolved determination of gaseous or volatile electrochemical reaction intermediates and products in real time.



A novel differential electrochemical mass spectrometry (DEMS) cell for integration with the Hidden Analytical HPR-40 DSA membrane inlet mass spectrometer systems. The new cell, developed by Ezra




► Figure 1. DEMS results obtained for CO<sub>2</sub>-sparged 0.05 M K<sub>2</sub>CO<sub>3</sub> electrolyte (pH = 6.8) with an electrolyte flow rate of 1 mL/min and a scan rate of 0.2 mV/s. Further details are included in the ACS publication.

L. Clark and Prof. Alexis T. Bell of The University of California, Berkeley Lab, coupled with Hidden's differentially pumped mass spectrometer system provides for in-situ mass resolved determination of gaseous or volatile electrochemical reaction intermediates and products in real-time.

## QIC series gas analysers with optimised sampling for real-time off-gas analysis

Hidden's QIC series gas analysers provide for multi-component, wide dynamic range real-time analysis of the key species involved in electrochemistry, hydrogen, oxygen, carbon dioxide and reaction products; ethanol for example.

A recent publication in Nature includes data from the Hidden MS: Jun Lu et al. (2016) "A lithium-oxygen battery based on lithium superoxide" Nature **378** (529), 377–382.



► Figure 2. Gas evolution results of Li-O<sub>2</sub> cells c) without a catalyst and d) with a catalyst while charging as measured by DEMS.




Figure 1: E. L. Clark, M. R. Singh, Y. Kwon, and A. T. Bell (2015) "Differential Electrochemical Mass Spectrometer Cell Design for Online Quantification of Products Produced during Electrochemical Reduction of CO<sub>2</sub>" Anal. Chem., **87** (15), 8013–8020

Figure 2: K. Kang et al. (2013) "Mechanism of Co<sub>3</sub>O<sub>4</sub>/graphene catalytic activity in Li-O<sub>2</sub> batteries using carbonate based electrolytes" Electrochimica Acta **90**, 63–70