


HPR-60

Molecular Beam Sampling Mass Spectrometer




- Reaction Kinetics
- Catalysis
- TPD/TPR/TPO
- TA-MS
- Gas Purity
- Environmental Gas Analysis
- Combustion
- CVD/MOCVD


HPR-60 detail



alignment bellows detail




application – carbon arc evaporation



mass spectrum from 300 - 1000 amu of Fullerene C₆₀ produced from high pressure arc evaporation. The zoomed series of C₆₀ peaks are identified as incorporating upto n=9 of the ¹³C isotope.


application – plasma sampling



schematic of plasma reactor and 3 stage differentially pumped line of sight mass spectrometer with cross beam ion source

Singh et al J. Vac. Sci. Tech. A17, 2447, (1999)


application – plasma sampling



beam to background measurements as a function of number of pumping stages and distance from skimmer to ion source

Singh et al J. Vac. Sci. Tech. A17, 2447, (1999)

application – plasma sampling



measurement
of plasma
radicals

Singh et al J. Vac. Sci. Tech. A17, 2447, (1999)

application – combustion studies




Fig. 1. Schematic diagram of the new flame-ion mass spectrometer.

flame
ionisation
mass
spectrometer

Goodings et al Int. J. Mass Spec. & Ion Proc. 132, 83, (1994)

www.HidenAnalytical.com

application – CVD diamond





FIG. 1. Schematic diagram of the molecular beam mass spectrometer coupled to the hot filament CVD reactor.

Rego et al Int. J. Appl. Phys.79, 7264, (1996)

www.HidenAnalytical.com

application – CVD diamond



stable gas species and methyl radicals as a function of filament temperature during CVD diamond deposition

Rego et al Int. J. Appl. Phys.79, 7264, (1996)

mass spectrometer options

	HAL/3F RC	HAL/3F PIC	EPIC	EQP
Detection Scheme	Analogue	Pulse Ion Counting	Pulse Ion Counting	Pulse Ion Counting
Measured	Neutrals Radicals	Neutrals Radicals	Neutrals Radicals + Ions - Ions retarding field energy analysis	Neutrals Radicals + Ions - Ions IED's

mass ranges to 2500 amu

ppb detection levels

MASsoft - [MID Tabular view C:\HIDEN\H2HPPMAR.EXP View:3]


File Edit MassSpecs Tune Gallery System Views SPC Window Help

Parameters/Scans:	ppm Hydrogen 1	ppm Air 1	ppm Water 1	ppm Methane 1	ppm Oxygen 1	ppm CO2 1
Real time						
7:57:32 pm	31	3	4	0.7	0.4	0.7
7:58:42 pm	32	2	3	0.5	0.04	0.6
7:59:55 pm	33	3	3	0.7	0.5	1
8:01:14 pm	34	2	3	0.5	0.5	0.9
8:02:27 pm	37	4	3	0.5	0.3	0.6
8:03:40 pm	36	2	4	0.5	0.08	0.8
8:04:53 pm	32	3	3	0.9	0.2	0.9
8:06:06 pm	35	2	3	0.4	0.3	0.8
8:07:19 pm	34	2	3	0.7	0.4	0.5
8:08:32 pm	35	1	5	0.8	0.4	0.7
8:09:57 pm	37	2	3	0.6	0.3	0.7
8:11:15 pm	33	3	2	0.7	0.2	0.8
8:12:28 pm	34	3	3	0.7	0.2	1
8:13:41 pm	34	3	4	0.8	0.2	0.6
8:15:00 pm	32	2	2	0.9	0.1	0.8
8:16:19 pm	32	2	4	0.6	0.2	0.5
8:17:38 pm	33	3	4	0.6	0.2	0.6
8:18:51 pm	36	2	3	0.8	0.5	0.7
8:20:04 pm	31	3	4	0.8	0.3	0.8
8:21:23 pm	32	2	3	0.8	0.2	0.8
8:22:36 pm	35	3	4	0.5	0.4	0.7
8:23:49 pm	33	3	4	0.8	0.3	1
8:25:02 pm	37	2	3	0.8	0.3	0.8
8:26:15 pm	33	2	4	0.7	0.2	0.7
8:27:40 pm	31	2	4	0.5	0.3	0.8
8:28:59 pm	33	2	4	0.4	0.3	0.8
8:30:12 pm	37	3	2	0.5	0.2	0.9
8:31:25 pm	29	2	2	0.8	0.3	0.7

compact HPR-60 MBMS




software examples



Scan Gallery

- simple 'Scan Gallery' provides rapid, easy to use set up of different types of scans.
- only 3 clicks of the mouse required to acquire mass spectra from MASsoft startup program.


software examples



Multiple Ion Detection Mode

- Scan up to 100 different masses simultaneously.
- Choose masses from the internal library.
- View data as graphical views, tabular views or both at the same time.
- Display output in preferred units (i.e. partial pressure, %, ppb)

software examples




Parameters/Scans:	ppm Krypton 1	ppm Helium 1	%Argon 1	ppm CO2 1
Mean	1.076	4.293	1.2	308.316
Range	1.163	3.014	0.0405	28.430
Standard Deviation	0.2546	0.5292	0.009191	5.198
Elapsed time				
16:38.476	0.6164	4.432	1.204	307.114
17:03	1.052	4.594	1.109	298.134
17:27.496	0.6865	4.813	1.211	310.071
17:52.050	0.9516	4.604	1.18	301.826
18:16.600	0.7315	3.958	1.196	308.912
18:41.100	1.01	4.71	1.197	308.99
19:05.532	0.8715	3.928	1.208	307.117
19:30.105	0.8821	4.430	1.108	302.595
19:54.585	1.591	4.711	1.19	308.155
20:19.886	0.9122	4.508	1.21	313.183
20:43.655	0.5779	4.152	1.196	304.950
21:08.190	1.236	3.671	1.195	305.331
21:32.669	1.018	3.33	1.194	306.393
21:57.294	1.155	4.349	1.19	311.249
22:21.734	0.8833	4.433	1.193	306.815
22:43.254	0.6761	3.94	1.185	305.698
23:07.831	1.366	4.624	1.203	305.325
23:32.350	1.227	3.534	1.214	303.014
23:56.887	1.349	3.554	1.201	315.646
24:21.424	1.266	3.835	1.212	311.106
24:45.982	1.341	4.668	1.192	314.469
25:10.481	1.386	5.443	1.202	306.451
25:34.924	1.082	3.745	1.208	304.476
25:59.425	0.991	4.298	1.188	302.531
26:23.972	0.7894	4.201	1.207	308.836

Tabular View

- results presented as a function of time in mixed units (ppm and %).
- the yellow markers indicate the range of data that has been selected for statistical analysis.
- Mean, Range and Standard Deviation from this analysis are shown at the top of the view and updated in real time.

software examples



I/O Events Editor

- powerful Events Editor allows trip, threshold and alarm levels to be set easily on screen. Output is to screen, audible beep or integral TTL and 0-10V analogue switches.
- used for 'out of limit' alarm indicators or sophisticated process control routines

software examples



LabVIEW™ drivers

- the Hiden drivers are automatically self-configuring, with all mass spectrometer control parameters grouped within LabVIEW™ for convenient programming for the full application range .

why choose Hiden Analytical?

company profile

Hiden Analytical was founded in 1981 and is presently situated in a 23,000 sq. ft. manufacturing plant in Warrington, England with a staff of 65 persons.

Hiden Analytical Inc, a wholly-owned subsidiary of Hiden Analytical Ltd, was formed in New Hampshire on January 1st 1996 to establish a domestic USA sales/service centre.



Certificate No. 6738

why choose Hiden Analytical?

people

- ◆ 65 staff providing...
- ◆ sales & service on 4 continents with ...
- ◆ 20 years manufacturing experience and...
- ◆ over 100 staff publications in peer reviewed journals
- ◆ over 200 user publications in peer reviewed journals

why choose Hiden Analytical?

Installations the following sites use Hiden Gas Analysis Systems

USA

Applied Materials

Axelis

CVC/Veeco

DuPont

General Motors

IBM Research

Lawrence Livermore

Motorola

NIST

Semtech

UK/Europe

Bosch

IMEC

Motorola

Nortel Networks

Oxford Plasma Technology

Philips

Rolls Royce

SGS Thomson

Siemens

Surface Technology Systems

Asia Pacific

Canon

Hitachi Fundamental Res.

Hyundai

LG Electronics

NEC

Samsung

Sony Corporation

TDK

Tokyo Electron

Toshiba