

Screening of catalysts for H₂S abatement from biogas to feed molten carbonate fuel cells

Our research project concerns the H₂S abatement from a biogas stream to feed to molten carbonate cells (MCFC) for energy production.

Biogas is a renewable energy source that can be produced from biomass, agricultural and industrial waste and sewage sludge. Its main constituents are CH₄, CO₂ but unfortunately sulphur based compounds are present as H₂S, COS.

The main limitation to the use of biogas as fuel in these systems is linked to the presence of sulphur compounds such as H₂S that can poison the fuel cells components, anode and electrolyte.

For the small scale plant, a very interesting solution can be represented by the direct and selective H₂S oxidation to sulphur at low temperature by using a suitable catalyst.

The aim of the work is the study of the catalytic behavior of V₂O₅/CeO₂ samples at various vanadium loading (2.55 - 20 wt % as V₂O₅) to determine the optimal catalyst formulation and operating conditions to obtain a very high H₂S conversion and the SO₂ selectivity minimisation.

The catalytic tests were carried out in a fixed bed flow reactor, made of a steel tube 21 cm long and a 14 mm of internal diameter. The reactor is inserted in an electrical furnace equipped with a PID electronic temperature controller. A thermocouple is inserted in a steel sheath of the inner diameter of 6 mm concentric to the reactor.

Catalytic activity tests were carried out at atmospheric pressure and GHSV=9.5·10⁴ h⁻¹ (40 ms), in the temperature range 150 – 250 °C, by feeding a mixture of 200 ppm of H₂S, 100 ppm of O₂ and N₂ to balance (see Fig. 1)




Fig. 1: Experimental apparatus for the catalytic tests

The exhaust system was equipped with a sulphur trap and analysed by a quadrupole mass spectrometer (Hidden Analytical, HPR-20 QIC). Figure 2 shows typical behaviour of a catalytic activity test carried out by using a catalyst using the highest vanadium loading (20%). With this condition being the optimal catalytic performance and more selective to sulphur (99%) at lower temperature.




Fig.2: Catalytic activity test of 20% V₂O₅/CeO₂ catalyst at 150 °C

It can be seen from the data above, in the first 10 minutes of the experiment, that the feed stream was added in the by-pass position to stabilise the initial concentrations of H₂S and O₂. After this time, the feed stream was sent to the reactor, and large increases in the concentrations of H₂S, SO₂ and O₂ were observed. After a further 10 minutes the concentrations stabilised, showing high activity and stability of the catalyst.

The H₂S and O₂ conversions were 97 and 98% respectively, while the SO₂ concentration remained very low (~1%).

Project summary by:

V. Palma & D. Barba
Department of Industrial Engineering
University of Salerno
Via Ponte Don Melillo
Fisciano (SA)
Italy



Paper Reference:

V. Palma, D. Barba, P. Ciambelli (2013) "Screening of catalysts for H₂S abatement from biogas to feed molten carbonate fuel cells" *International Journal of Hydrogen Energy* **38** (1) 328-335

Hidden Product:

HPR-20 QIC Real time Gas Analyser

Follow the link to the product catalogue on our website for further information

<http://www.hiddenanalytical.com/en/products/for-gas-analysis/hpr-20-qic-r-d#customer-stories>