Differentially pumped quadrupole SIMS probe on FIBbased and two-beam microscopes




UK


Richard J Chater⁽¹⁾, Barbara Shollock⁽¹⁾, David McPhail⁽¹⁾, Alan J Smith⁽²⁾ and Graham Cooke⁽²⁾

• (1) Department of Materials, Imperial College London, SW7 2AZ, UK


Imperial College London

• (2) Hiden Analytical Ltd, Warrington, WA5 7UN,









Zeiss Auriga SEM

 FEG SEM electron gun column with electrostatic final lens.

♦ Gallium ion beam column
 from Orsay-Physics (Cobra) with
 beam spot sizes that vary from
 10nm.

- Secondary electron detectors both in-lens and separate
 - Backscatter electron detectors in-lens
- Positive or negative charged particle detector (SESI)

♦ EDS

♦ EBSD

Hiden EQS 1000 SIMS detector

 ◆ Triple quadrupole electric mass filter for masses from 0.4 AMU to 300AMU

- Electrostatic filter for ions at quadrupole entrance.
- secondary charged ions detected individually by secondary electron multiplier.
- separate vacuum pumping using a drypump and turbomolecular pump.
- Software system for detector setup and control for spectra, depth profiles and images.
- Residual gas analysis (RGA)


FEI FIB200 SIMS


 FIB200 workstation with single beam gallium ion gun used at energies to 30keV.
 Beam can be scanned with normal line/frame raster or within a pattern(s).


 ♦ Gallium ion beam spot size varies from 10nm at 20pA to ~600nm at 20nA.

◆ FEI designed quadrupole based SIMS detector with low field collection and without an electrostatic analyser in the secondary ion column.

◆ SIMS measurements generate spectra, depth profiles or maps.







Vibration effects from image analysis

♦ Gold clusters on carbon
 imaged with Gallium beam at
 30keV and 100pA aperture in
 both the FEI FIB200 and Zeiss
 Auriga SEM


◆ Line intensity profiles horizontally and vertical through the centre of the images were analysed for the power spectra.

 Drypump and turbopump damping and connections modified in response to power spectra

FEI FIB200


 Image degradation is not seen at any beam current after the modification to the Hiden SIMS probe pumping system.


◆ Zeiss Auriga SEM-FIB





 Vibrational effects of the additional pumping system with the Hiden SIMS probe is seen in the SEM image of the gold cluster on carbon sample at 200kX magnification.

 Further damping for low frequency is required









Pressure differential : FEI FIB200

- ◆ FEI FIB200 instrument has a system for introducing water
 vapour only for positive ion yield M⁺ enhancement.
- Chamber pressure is controlled by the reservoir temperature which contains hydrated magnesium sulphate (Epsom Salts).
- Graph shows that a differential pressure of ~2 orders of magnitude can be maintained.
- Hiden SIMS probe would have tripped off at the lowest reservoir temperature of 25°C.

Pressure differential : Zeiss Auriga SEM

- Zeiss Auriga SEM has a system for introducing oxygen gas as a jet directed at the sample surface for charge compensation. SE and BSE ionise local gas molecules to achieve a charge balance for normal imaging.
- Chamber pressure is controlled by flowrate and gas source pressure.
- Approximately two orders of magnitude pressure differential is maintained.

- Yield enhancement ratio for M⁺ with oxygen coverage compared to a clean surface.
 - ◆ A. Benninghoven 1976

Z	М	R
24	Chromium	1000
38	Strontium	800
42	Molydenum	615
25	Manganese	500
74	Tungsten	389
22	Titanium	308
23	Vanadium	300
73	Tantalum	286
26	Iron	233
56	Barium	150
13	Aluminium	100
12	Magnesium	90
41	Niobium	83
28	Nickel	75
14	Silicon	69
29	Copper	23
32	Gemanium	5

Platinum layer sectioned at crater wall and visible by tilting the sample


 PÅ
 Mag
 Det
 Scan
 HFW
 FWD
 Tit
 1 µm

Depth profile of Boron in Silicon


- Small area, 70um x 70 um coated in-situ with ~200nm of a platinum organic.
- Ga⁺ FIB sputtering (30keV, 3nA) into area of 50um x 50 um centered within the platinum.
- Inset image shows the crater wall and crater base for depth estimation.
- Water vapour in main chamber, pressure 8.8 x 10⁻⁶ mbar.
- Platinum remains at edge throughout the depth profiling despite gallium beam skirt.

¹¹B implant, dose 1.6 x 10¹⁵ per cm² at 40keV into silicon

- Peak concentration ~0.25%at
- SRIM estimates ¹¹B peak at 144nm depth below surface which is distinguishable on the Log-Lin plot.
- Oxygen enhancement for both silicon and boron are below 100.
- TofSIMS profile, Oxygen sputtering at 1keV. 25keV Bi pulses.



TofSIMS measurement of Depth Profile




Simultaneous detection of positive and negative secondary ions in the single beam FEI FIB200 SIMS instrument


Lanthanum strontium Manganate target

Positive SIMS_FEI SIMS detector

Negative SIMS_HIDEN SIMS detector

- Gallium FIB ion beam at 30keV and 20nA into a crater of 50 microns square on the surface of the target.
- Chamber pressure during sputtering was 7.3 x 10⁻⁷ mbar using a water vapour leak into the chamber for positive ion enhancement.
- Two separate mass spectra results are shown in the chart with simultaneous SIMS detection whilst the single crater was sputtered into the target.





Particle emission for low energy ion impacts

 Positive and negative secondary ion detection using both the Hiden SIMS probe and the Zeiss SESI detector.

 ♦ Secondary electron detection using inlens secondary electron detectors within the Zeiss electron column and the Zeiss SESI detector.

 ♦ X-Ray emission detector using standard Oxford Instruments SDD EDS/EDX detector.

Discussion

- Feature selection is often best achieved by exploiting a range of well established and very familiar analytical techniques that are standard on SEM instruments.
- Correct sample mounting and orientation in different instruments from different manufacturers can be very difficult. Finding micron-sized features on transfer between instruments can be very time consuming.
- Recent developments in low energy electron beam columns now allow for both electron and ion beams to have approximately the same excitation volume on the sample surface by correct selection of electron energy and ion energy. This leads to imaged data at nanometer resolution that is directly comparable. SIMS compositional information can be matched to the topographical and structural information from secondary electron and EBSD imaging.

Conclusions

- The added SIMS analytical facility to both the two beam and single beam microscopes has demonstrated its potential in a very short time since its installation in late July 2013.
- SIMS facilities alongside the more mainstream analytical techniques available on most two beam microscopes ensures an increased access and awareness of the technique.
- Potential for simultaneous positive and negative SIMS has been achieved.

Acknowledgements

Mahmood Ardakani and Dani Proprentner of Imperial College and the splendid team at Hiden Analytical for their enthusiam for this project also Dr Giles Graham , AWE, UK for his encouragement and vision .