

## **FIB-SIMS**

Powerful Surface Analysis with Focused Ion Beam -  
Secondary Ion Mass Spectrometry Depth Profiling



# FIB-SIMS

## Powerful Surface Analysis with Focused Ion Beam - Secondary Ion Mass Spectrometry Depth Profiling

---

Contact: Dr. Ulrich Kohl-Roscher

---

Date: March 2016

---

**Analyze the composition of solid surfaces and thin films by sputtering the surface of the specimen.**

### Introduction

FIB-SIMS is a very powerful surface analytical technique especially for high sensitivity nano-scale materials analysis. Elemental detection limits range from parts per million to parts per billion. Generate elemental surface, image and depth profile information by mass spectrometry. Use the FIB-SIMS to analyze the composition of solid surfaces and thin films with a focused primary ion beam and collecting and analyzing ejected secondary ions. SIMS is considered to be a qualitative technique, although quantitation is possible with the use of standards.

### Availability

The FIB-SIMS is available for the following microscopes:

- AURIGA series
- NEON series
- 15xx series
- XB 340, XB 540


### Benefits

- Analysis of trace elements down to ppm (parts per million) level (thin films, semiconductors, solar cells)
- Isotope detection
- Elemental mapping and depth profiling
- Detection of atomic and molecular ions
- 50 nm lateral analytical resolution possible (better than conventional SIMS)
- A primary beam of energetic ions, typically in the range of 500 V up to 30 kV is used to sputter or erode the surface of the material under analysis

### Operation

A primary focused ion beam is interacting with the surface and generates charged (electron, ions) and uncharged (molecules and elements) particles (see fig. 1). The sputtered charged ions / molecules are directed into the detector by electric fields. The ions pass through an energy analyzer and through a mass analyzer. The mass and energy resolved ions are detected using a secondary electron multiplier. Sputtered ions from milling are collected and analyzed by mass spectroscopy. The detected ions provide analytical information. SIMS is a very surface sensitive technique because secondary ions emit only from the upper monolayers of the sample.

Figure 1



**Upgrade path**

Software	All SmartSEM versions
Hardware	SIMS port required
	FIB column

A technical feasibility check has to be performed before issuing a quotation.

A system preventive maintenance performed within the last 12 months is mandatory.

The retrofit must be performed by a ZEISS-authorized service engineer.

Application training is recommended.

For further information, contact:

[microscopy@zeiss.com](mailto:microscopy@zeiss.com)

Part	Ordering no.
Upgrade Kit SIMS EQS	346565-8002-990
Upgrade Kit SIMS CUST XB	346565-8003-990
Upgrade Kit SIMS CUST AC	346565-8004-990
Installation and on site training (3 day system installation and training at customer site provided by Hiden Analytical Ltd.), Europe	346565-8005-000
Installation and on site training (3 day system installation and training at customer site provided by Hiden Analytical Ltd.), USA & Canada	346565-8006-000
Installation and on site training (3 day system installation and training at customer site provided by Hiden Analytical Ltd.), Far East	346565-8007-000
Additional 12 months warranty extension (cumulative up to a maximum of 5 years total)	346565-8008-000



**Carl Zeiss Microscopy GmbH**

07745 Jena, Germany  
microscopy@zeiss.com  
www.zeiss.com/microscopy



We make it visible.