


## Unexpected Low-Temperature DeNO<sub>x</sub> Activity of AdSCR Systems for Cold start NO<sub>x</sub> Abatement

Urea-SCR aftertreatment systems have been implemented in lean burn vehicles to reduce NO<sub>x</sub> emissions in order to meet government regulations. Current commercial catalysts (Cu- and Fe- promoted zeolites) are able to achieve very high NO<sub>x</sub> removal efficiencies once they reach operating temperatures above 200°C. However, their performances are relatively poor at low temperatures, e.g. during the engine cold start. Indeed, in this temperature region: (i) the inability of decomposing urea prevents fully exploiting the SCR technology; (ii) the metal-promoted zeolite catalysts have a limited deNO<sub>x</sub> activity.

In the present work we extend the investigation of a novel catalytic device designed to improve the abatement of NO<sub>x</sub> during the engine warm up phase, the so called AdSCR system (AdSCR = Adsorption + Selective Catalytic Reduction). The formulation of AdSCR systems is based on the combination of a Urea-SCR catalyst (e.g. Cu-CHA) and of a NO<sub>x</sub> storage material (e.g. BaO/Al<sub>2</sub>O<sub>3</sub>) in the form of physical mixtures. They are characterized by a dual functionality: they can operate both as NO<sub>x</sub> adsorbers, trapping NO<sub>x</sub> during the low temperature transient, and as SCR catalysts, reducing the previously stored NO<sub>x</sub> with NH<sub>3</sub> at higher temperatures.

A set of experiments was carried out over a physical mixture of SCR catalysts (Cu-CHA) and NO<sub>x</sub> storage materials (BaO/Al<sub>2</sub>O<sub>3</sub>) to study their deNO<sub>x</sub> efficiency. The typical protocol consisted of cold start mimicking runs, aimed at simulating the aftertreatment system warm-up: NO (300 ppm) was fed in a stream of 8% of O<sub>2</sub> at r.t. followed by a temperature (T)-ramp (15 °C/min). At 170 °C, i.e. the T-threshold for urea injection, NH<sub>3</sub> (800 ppm) was fed to start the SCR reactions. The same protocol was repeated in the presence of NH<sub>3</sub> pre-adsorbed at 320 °C to investigate the reactivity of stored NO<sub>x</sub> and NH<sub>3</sub>. In these experiments the dynamics of NO, NH<sub>3</sub>, and NO<sub>2</sub> at the reactor outlet were continuously measured both by a UV analyzer (ABB LIMAS 11HW) and by a quadrupole mass spectrometer (QGA Hidden Analytical) which could, furthermore, detect N<sub>2</sub> and N<sub>2</sub>O.

Our results from dry cold-start runs show that: i) after a NO step feed, NO<sub>x</sub> were completely trapped for up to 420 s before NO breakthrough, with a NO storage efficiency of 93% before NH<sub>3</sub> injection; ii) after NH<sub>3</sub> injection, a drop in the NO outlet concentration and a release of N<sub>2</sub> (850 ppm) were observed due to the onset of the SCR reactions. Moreover, in the presence of pre-adsorbed ammonia the formation of N<sub>2</sub> was unexpectedly observed already at 40 °C. These remarkable results identify AdSCR as promising systems and may warrant further investigation.



*The Hiden QGA in the Laboratory*

### Project summary by:

Enrico Tronconi, Laboratory of Catalysis and Catalytic Processes, Dipartimento di Energia  
Politecnico di Milano, Via La Masa 34  
20156 Milano, Italy



**POLITECNICO**  
MILANO 1863

### Paper Reference:

Federica Gramigni, Nicole Daniela Nasello, Tommaso Selleri, Isabella Nova, Enrico Tronconi, Simone Dieterich, Michel Weibel "Unexpected Low-Temperature DeNO<sub>x</sub> Activity of AdSCR Systems for Cold start NO<sub>x</sub> Abatement" *Emiss. Control Sci. Technol.* (2020) 6:402–409

### Hidden Product:

[QGA](#)