


A Robust Two-dimensional Layered Metal-Organic Framework for Efficient Separation of Methane from Nitrogen

Fang Zheng^a, Lihang Chen^{a, *}, Rundao Chen^a, Zhiguo Zhang^{a, b}, Qiwei Yang^{a, b}, Yiwen Yang^{a, b}, Baogen Su^{a, b}, Qilong Ren^{a, b}, Zongbi Bao^{a, b, *}

^aKey Laboratory of Biomass Chemical Engineering of the Ministry of Education, College of Chemical and Biological Engineering, Zhejiang University, 38 Zheda Road, Hangzhou 310027, P. R. China.

^bInstitute of Zhejiang University-Quzhou, Quzhou, 78 Jiu Hua Boulevard North, Quzhou 32400, P. R. China.

Selective capture and separation of methane (CH₄) from nitrogen (N₂) is a feasible approach to mitigate the effects of global warming and to improve the heating value of low-quality natural gas. However, the efficient separation of CH₄/N₂ is a challenge issue since their very close kinetic diameters and thermodynamic properties. In this work, we successfully combine twisted organic ligand (4-DPDS = 4,4'-dipyridyldisulfide) and angular inorganic anion (CrO₄²⁻) forming a new class of robust two-dimensional (2D) layered rhomboid grids with ionic polar environment. This a new type of 2D layered metal-organic framework (MOF) with permanent rhomboid pore channels, denoted as Ni(4-DPDS)₂CrO₄, was synthesized for the first time and used for CH₄/N₂ separation. Such a 2D layered material is rarely used for methane and nitrogen separation and shows high affinity towards CH₄ thereby affording admirable methane adsorption and separation performance. Specifically, it not only performs high stability as outstanding as or even better than previously reported 3D MOFs but also exhibits a relatively high CH₄ adsorption capacity of 0.95 mmol/g. Ni(4-DPDS)₂CrO₄ has a high affinity towards CH₄ with the highest reported Q_{st} value of 28.4 kJ/mol as well as a high CH₄/N₂ selectivity of 7.3, which is comparable to the state-of-the-art MOF materials reported so far.



Two-dimensional layered MOF

By DFT-D calculations, it was found that twisted organic ligand plays a vital role in delivering a relatively rigid structure and the angular inorganic anions provides polar sites for strong guest-host interactions in close proximity. It revealed the energy favorable binding sites for methane molecules are located in the middle of the cavity decorated with CrO₄²⁻ anion, where inducing C-H...O hydrogen bond with the distance of 2.52 Å. Simultaneously, each adsorbed CH₄ molecule is surrounded by eight pyridyl rings, among which one forms dominant cooperative supramolecular interactions

between C(δ^-) of CH₄ and H(δ^+) from the pyridyl ring (C \cdots H = 3.43 Å) and moderate-to-weak $\sigma\cdots\pi$ stacking interactions between the electrons in the π orbitals of the pyridyl ring and the electrons in H of CH₄ (H \cdots pyridyl ring = 3.47 Å). Furthermore, the weak intermolecular electrostatic dipole interactions between the H(δ^+) of CH₄ and S (δ^-) (H \cdots S = 3.82 Å) are also observed. Such an angular inorganic anion provides a cooperative effect conferred by polar sites and multiple guest-host interactions in close proximity, affording tight binding affinity.

Dynamic breakthrough experiments were carried in a stainless-steel column manually packed with Ni(4-DPDS)₂CrO₄. The outlet gas passing through the column was analyzed using a **Hidden HPR-20 EGA mass spectrometer** for continuous sampling gas analysis. It certified N₂ was eluted almost immediately and then quickly approached a pure grade while no detectable CH₄, whereas CH₄ was retained in the packed column for a remarkable time (~15 min) until reaching the saturated uptake and breakthrough. After the concentration of eluting gas remained unchanged, the amount of CH₄ adsorbed into the column reached 0.94 mmol/g, with a selectivity of 4.7 to N₂ (0.2 mmol/g), which suggest it is a promising material in natural gas purification by further regenerability tests. It is a thriving exploration in two-dimensional MOFs that will emerge the development of various materials for CH₄ trap in the future.

Project summary by:

Zongbi Bao

Key Laboratory of Biomass Chemical Engineering of the Ministry of Education,
College of Chemical and Biological Engineering,
Zhejiang University,
38 Zheda Road,
Hangzhou 310027
P. R. China.



Paper Reference:

"A robust two-dimensional layered metal-organic framework for efficient separation of methane from nitrogen" *Separation and Purification Technology* (2022) 281, 119911

Hidden Product:

[HPR-20 EGA](#)