


## Control of $Ti_{1-x}Si_xN$ nanostructure via tunable metal-ion momentum transfer during HIPIMS/DCMS co-deposition

Low-energy inert-gas ion irradiation of the film surface during refractory transition-metal (TM) nitride growth by conventional DC magnetron sputtering has been used extensively to overcome the characteristically underdense microstructures with rough surfaces of layers deposited at low temperatures ( $T_s/T_m < 0.30$ , in which  $T_s$  is the film growth temperature and  $T_m$  is the melting point in K).<sup>1</sup> We recently demonstrated that high-power pulsed magnetron sputtering (HIPIMS) provides an alternative route for ion-assisted TM nitride film growth via the use of substrate bias, synchronized to the *metal-rich* portion of the plasma pulse. Stresses can be dramatically reduced, or even eliminated, since *metal* (as opposed to inert-gas) ions are components of the film.<sup>ii,iii</sup>

In this project we use a hybrid HIPIMS/DCMS two-target co-sputtering configuration, in which one target (either Ti or Si) is powered by HIPIMS while the other is powered by DCMS, for growth of  $Ti_{1-x}Si_xN$  films with compositions  $0 \leq x \leq 0.26$ . Markedly different film growth pathways are obtained depending upon which target is powered by HIPIMS with, in both cases, a substrate bias applied in synchronous with the HIPIMS pulse. The observed divergence in film nanostructure, phase content, and mechanical properties between layers grown in Ti-HIPIMS/Si-DCMS and Si-HIPIMS/Ti-DCMS configuration is due to distinctly different metal-ion irradiation conditions,  $Ti^+/Ti^{2+}$  vs.  $Si^+/Si^{2+}$ , during film growth, as determined by the ion mass spectrometry analyses performed at the substrate position with a **Hidden Analytical EQP 1000 instrument** (see Fig. 1(a)-(b)).



**Figure 1.** Ion energy distribution functions measured at the substrate position for (a) singly-charged  $Ti^+$  and  $Si^+$  ions, and (b) doubly-charged  $Ti^{2+}$  and  $Si^{2+}$  ions during Ti-HIPIMS and Si-HIPIMS pulses; (c) plan-view STEM micrograph, and (d) plan-view EDX/STEM elemental maps of a  $Ti_{0.74}Si_{0.26}N$  Ti-HIPIMS/Si-DCMS film, showing spatial distributions, acquired from the area outlined in panel (c); (e) nanoindentation hardnesses  $H(x)$  of Ti-HIPIMS/Si-DCMS and Si-HIPIMS/Ti-DCMS  $Ti_{1-x}Si_xN$  films grown on Si(001) substrates at  $T_s = 500^\circ C$ .

A better mass match between incident  $Ti^+$  ions and the average film atomic mass, higher metal-ion/metal-atom ratios, and a high fraction of doubly-ionized species results in an average momentum transfer per deposited atom ( $\langle p_d \rangle$ ) ~20 times higher for Ti-HIPIMS/Si-DCMS than during Si-HIPIMS/Ti-DCMS. As a consequence, adatom mean free paths are increased leading to the segregation of smaller Si atoms to column boundaries and the formation of a nanocomposite structure consisting of TiN-rich nanocolumns encapsulated in  $SiN_x$  matrix phases (cf. plan-view STEM micrograph in Fig.1(c), and EDX/STEM elemental maps in Fig. 1(d)). Ti-HIPIMS/Si-DCMS  $Ti_{1-x}Si_xN$  films are superhard over a composition range that is significantly wider than reported previously,  $0.04 \leq x \leq 0.26$ , with a maximum hardness,  $H = 45$  GPa, for layers with  $x = 0.13$  (see Fig. 1(e)). However, residual stresses are also high with an average value of  $-7 \pm 1$  GPa.

In sharp contrast, during Si-HIPIMS/Ti-DCMS  $Ti_{1-x}Si_xN$  film growth, the flux of doubly-ionized metal ions is lower which, together with the lower mass of Si, low metal-ion/metal-atom flux ratio during HIPIMS pulses, and poorer mass match between incident  $Si^+$  ions and average film atomic mass results in relatively low ( $\langle p_d \rangle$ ) values. As a consequence, Si is trapped in the metastable  $Ti_{1-x}Si_xN$  NaCl structure to form solid solutions over the highest compositional range yet reported,  $0 \leq x \leq 0.24$ .

#### Project summary by:

G. Greczynski  
Thin Film Physics Division  
Department of Physics (IFM)  
Linköping University  
SE-581 83 Linköping  
Sweden



#### Paper Reference:

G. Greczynski et al., (2015) "Control of  $Ti_{1-x}Si_xN$  nanostructure via tunable metal-ion momentum transfer during HIPIMS/DCMS co-deposition", *Surface and Coatings Technology* **280**, 174-184  
<https://doi.org/10.1016/j.surfcoat.2015.09.001>

#### Hidden Product:

[EQP](#)

---

<sup>i</sup> I. Petrov, P.B. Barna, L. Hultman, J.E. Greene *J. Vac. Sci. Technol. A* 21 (2003) 117

<sup>ii</sup> G. Greczynski, J. Lu, I. Petrov, J.E. Greene, S. Bolz, W. Kölker, Ch. Schiffers, O. Lemmer and L. Hultman, *J. Vac. Sci. Technol. A* 32 (2014) 041515.

<sup>iii</sup> G. Greczynski, J. Lu, J. Jensen, I. Petrov, J.E. Greene, S. Bolz, W. Kölker, Ch. Schiffers, O. Lemmer and L. Hultman, *J. Vac. Sci. Technol. A* 30 (2012) 061504.